



TEST TUBES IN POLYSTYRENE WITH RIM

Tubos de ensayo en poliestireno con borde / Tubes en polystyrène avec bord / Provette in polistirolo con bordo



GRADUATED TEST TUBE
tubo graduado / tube gradué /
provetta graduata

MAXIMUM TRANSPARENCY
máxima transparencia / transparence maximale /
massima trasparenza

EXCELLENT RESISTANCE
excelente resistencia / excellente résistance /
eccellente resistenza



POLYSTYRENE
poliestireno / polystyrène / polistirolo

	CODE	CAPACITY	TYPE	PACKAGING
	21057	5 ml	13x75 mm (cylindrical)	IN BULK a granel / en vrac / sfusa
	21053	5 ml	12x86 mm (cylindrical)	IN BULK a granel / en vrac / sfusa
	21055	10 ml	16x100 mm (cylindrical)	IN BULK a granel / en vrac / sfusa
	21056	10 ml	16x100 mm (conical)	IN BULK a granel / en vrac / sfusa



TEST TUBES IN POLYSTYRENE WITHOUT RIM

Tubos de ensayo en poliestireno sin borde / Tubes en polystyrène sans bord / Provette in polistirolo senza bordo



HIGHER THICKNESS TEST TUBE
tubo de ensayo con espesor aumentado /
tube avec épaisseur augmenté /
provetta con spessore maggiorato

MAXIMUM TRANSPARENCY
máxima transparencia / transparence maximale /
massima trasparenza

EXCELLENT RESISTANCE
excelente resistencia / excellente résistance /
eccellente resistenza



POLYSTYRENE
poliestireno / polystyrène / polistirolo

	CODE	CAPACITY	TYPE	PACKAGING
	21085	10 ml	16x100 mm (cylindrical)	IN BULK a granel / en vrac / sfusa
	21086	10 ml	16x100 mm (conical)	IN BULK a granel / en vrac / sfusa