



SODIUM CITRATE TUBES FOR ESR

Tubos con Citrato de Sodio para VSG / Tubes avec Citrate de Sodium pour VS /
Provette con Sodio Citrato per VES



SEDI-RATE TEST TUBE WITH 0,2 ML OF SODIUM CITRATE (3,8%) X 0,8 ML OF BLOOD (ESR)

tubo Sedi-Rate con 0,2 ml de Citrato de Sodio (3,8%) x 0,8 ml de sangre (VSG) / Sedi-Rate tube avec 0,2 ml de Citrate de Sodium (3,8 %) x 0,8 ml de sang (VS) / provetta Sedi-Rate con 0,2 ml di Sodio Citrato (3,8%) x 0,8 ml di sangue (VES)

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite /
tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile

PINK

rosa / rose / rosa



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène /
polipropilene



TEST TUBE WITH 0,25 ML OF SODIUM CITRATE (3,8%) X 1 ML OF BLOOD (ESR) AND 2,25 ML OF BLOOD (COAGULATION)

tubo con 0,25 ml de Citrato de Sodio (3,8%) x 1 ml de sangre (VSG) y 2,25 ml de sangre (Coagulación) / tube avec 0,25 ml de Citrate de Sodium (3,8 %) x 1 ml de sang (VS) y 2,25 ml de sang (Coagulation) / provetta con 0,25 ml di Sodio Citrato (3,8%) x 1 ml di sangue (VES) e 2,25 ml di sangue (Coagulazione)

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite /
tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile

PINK

rosa / rose / rosa



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène /
polipropilene

| CODE | CAPACITY | TYPE | PACKAGING |
|-------|----------|----------|---|
| 50335 | 0,8 ml | 12x86 mm | RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi |

| CODE | CAPACITY | TYPE | PACKAGING |
|-------|--|----------|---|
| 22213 | 1 ml x ESR 2,25 ml x COAGULATION | 12x56 mm | RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi |
| 22218 | 1 ml x ESR 2,25 ml x COAGULATION | 13x75 mm | RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi |
| 22210 | 1 ml x ESR 2,25 ml x COAGULATION | 12x86 mm | RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi |