




SODIUM CITRATE TUBES FOR ESR

Tubos con Citrato de Sodio para VSG / Tubes avec Citrate de Sodium pour VS / Provette con Sodio Citrato per VES



TEST TUBE WITH 0,4 ML OF SODIUM CITRATE (3,8%) X 1,6 ML OF BLOOD (ESR)

tubo con 0,4 ml de Citrato de Sodio (3,8%) x 1,6 ml de sangre (VSG) / tube avec 0,4 ml de Citrate de Sodium (3,8 %) x 1,6 ml de sang (VS) / provetta con 0,4 ml di Sodio Citrato (3,8%) x 1,6 ml di sangue (VES)

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile



POLYETHYLENE


polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
--	------	----------	------	-----------




22212

1,6 ml x ESR

12x56 mm
(pink)

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



22257

1,6 ml x ESR

13x75 mm
(pink)

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



22209

1,6 ml x ESR

12x86 mm
(pink)

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



32267

1,6 ml x ESR

13x75 mm
(black)

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



SODIUM CITRATE TUBES FOR COAGULATION

Tubos con Citrato de Sodio para Coagulación / Tubes avec Citrate de Sodium pour Coagulation / Provette con Sodio Citrato per Coagulazione



TEST TUBE WITH SODIUM CITRATE FOR COAGULATION

tubo con Citrato de Sodio para Coagulación / tube avec Citrate de Sodium pour Coagulation / provetta con Sodio Citrato per Coagulazione

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile

YELLOW

amarillo / jaune / giallo



POLYETHYLENE


polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
--	------	----------	------	-----------




22259

0,25 ml of Sodium x
2,25 ml of blood

12x56 mm

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



22248

0,25 ml of Sodium x
2,25 ml of blood

13x75 mm

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



22250

0,4 ml of Sodium x
3,6 ml of blood

13x75 mm

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



22260

0,25 ml of Sodium x
2,25 ml of blood

12x86 mm

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



22255

0,4 ml of Sodium x
3,6 ml of blood

12x86 mm

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi



22254

0,5 ml of Sodium x
4,5 ml of blood

12x86 mm

RACK OF 50 PIECES
Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi