



PAEDIATRIC TUBES

Tubos pediátricos / Tubes pédiatriques / Provette per prelievo pediatrico



TEST TUBE WITH LITHIUM HEPARIN X 1 ML OF BLOOD
 tubo con Heparina de Litio x 1 ml de sangre /
 tube avec Héparine de Lithium x 1 ml de sang /
 provetta con Litio Eparina x 1 ml di sangue

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite /
 tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène /
 polipropilene



TEST TUBE WITH 0,1 ML OF SODIUM CITRATE (3,8%) X 0,9 ML OF BLOOD (COAGULATION)

tubo con 0,1 ml de Citrato de Sodio (3,8%) x 0,9 ml de
 sangre (Coagulación) / tube avec 0,1 ml de Citrate
 de Sodium (3,8%) x 0,9 ml de sang (Coagulation) /
 provetta con 0,1 ml di Sodio Citrato (3,8%) x 0,9 ml di
 sangue (Coagulazione)

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite /
 tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile

YELLOW

amarillo / jaune / giallo





POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène /
 polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
	22307	1 ml	12x56 mm (blue)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32307	1 ml	12x56 mm (dark green)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi

CODE	CAPACITY	TYPE	PACKAGING
22239	0,9 ml	12x56 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi