



SERUM SEPARATOR TUBES

Tubos separadores de suero / Tubes avec séparateur de sérum / Provette con separatore di siero





TEST TUBE WITH GEL SEPARATOR AND CLOT ACTIVATOR SILICONE COATED INTERIOR
 tubo con gel separador y activador de la coagulación interior recubierto de silicona / tube avec gel séparateur et activateur de la coagulation silicone à l'intérieur / provetta con gel separatore e acceleratore della coagulazione siliconata internamente

PERFECT TRACEABILITY
 trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

UNBREAKABLE
 irrompible / incassable / infrangibile

RED
 rojo / rouge / rosso

 **POLYETHYLENE**
 polietileno / polyéthylène / polietilene


 **POLYPROPYLENE**
 polipropileno / polypropylène / polipropilene





TEST TUBE WITH GEL SEPARATOR AND CLOT ACTIVATOR
 tubo con gel separador y activador de la coagulación / tube avec gel séparateur et activateur de la coagulation / provetta con gel separatore e acceleratore della coagulazione

PERFECT TRACEABILITY
 trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

UNBREAKABLE
 irrompible / incassable / infrangibile

 **POLYETHYLENE**
 polietileno / polyéthylène / polietilene

 **PMMA**
 metacrilato / méthacrylate / metacrilato

	CODE	TYPE	PACKAGING
	22744	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22734	12x86 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22702	16x100 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi

	CODE	TYPE	PACKAGING
	22742	13x75 mm (red)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22732	12x86 mm (red)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22722	16x100 mm (red)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22842	13x75 mm (brown)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22832	12x86 mm (brown)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	22822	16x100 mm (brown)	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi