




URINE CONTAINERS 120 ML


Frascos 120 ml para orina / Pots 120 ml pour urine / Contenitori urina 120 ml



URINE CONTAINER WITH FROSTED LABEL
 frasco para orina con banda mate / pot pour urine avec surface d'écriture / contenitore per urina con superficie di scrittura

WITH SCREW CAP APART
 con tapón no roscado / avec bouchon à vis no jointé / con tappo a vite a parte

 **POLYETHYLENE**
 polietileno / polyéthylène / polietilene


 **POLYPROPYLENE**
 polipropileno / polypropylène / polipropilene




URINE CONTAINER WITH FROSTED LABEL
 frasco para orina con banda mate / pot pour urine avec surface d'écriture / contenitore per urina con superficie di scrittura

WITH SCREW CAP
 con tapón roscado / avec bouchon à vis jointé / con tappo a vite

 **POLYETHYLENE**
 polietileno / polyéthylène / polietilene

 **POLYPROPYLENE**
 polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	STERILITY	PACKAGING
	25030	120 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25031	120 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa

	CODE	CAPACITY	STERILITY	PACKAGING
	25033	120 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25034	120 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25036	120 ml	STERILE estéril / stérile / sterile	IN BULK a granel / en vrac / sfusa
	25032	120 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	INDIVIDUALLY WRAPPED envuelto individualmente / sachet unitaire / in sacchetto singolo
	25035	120 ml	STERILE estéril / stérile / sterile	INDIVIDUALLY WRAPPED envuelto individualmente / sachet unitaire / in sacchetto singolo