



URINE CONTAINERS 60 ML

Frascos 60 ml para orina / Pots 60 ml pour urine / Contenitori urina 60 ml



URINE CONTAINER WITH FROSTED LABEL

frasco para orina con banda mate / pot pour urine avec surface d'écriture / contenitore per urina con superficie di scrittura

WITH SCREW CAP APART

con tapón no roscado / avec bouchon à vis no jointé / con tappo a vite a parte



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene



URINE CONTAINER WITH FROSTED LABEL

frasco para orina con banda mate / pot pour urine avec surface d'écriture / contenitore per urina con superficie di scrittura

WITH SCREW CAP

con tapón roscado / avec bouchon à vis jointé / con tappo a vite



POLYETHYLENE






polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	STERILITY	PACKAGING
	25180	60 ml (white)	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25181	60 ml (red)	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa

	CODE	CAPACITY	STERILITY	PACKAGING
	25183	60 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25184	60 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25186	60 ml	STERILE estéril / stérile / sterile	IN BULK a granel / en vrac / sfusa
	25182	60 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	INDIVIDUALLY WRAPPED envuelto individualmente / sachet unitaire / in sacchetto singolo
	25185	60 ml	STERILE estéril / stérile / sterile	INDIVIDUALLY WRAPPED envuelto individualmente / sachet unitaire / in sacchetto singolo