

CONTAINERS IN SINGLE BOX

Frascos en estuche / Pots en étui / Linea farmacia



URINE CONTAINER

frasco para orina / flacon carré pour urine / contenitore per raccolta urina

STERILE

estéril / stérile / sterile



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

CODE	CAPACITY	PACKAGING
25177	120 ml	INDIVIDUALLY WRAPPED, IN SINGLE BOX en bolsa individual, en estuche / emballage individuel, en étui / in sacchetto singolo, in astuccio

FAECES CONTAINER

frasco para heces / pot pour selles / contenitore per feci

STERILE

estéril / stérile / sterile



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

CODE	CAPACITY	PACKAGING
25193	60 ml	INDIVIDUALLY WRAPPED, IN SINGLE BOX en bolsa individual, en estuche / emballage individuel, en étui / in sacchetto singolo, in astuccio



URINE TEST TUBE WITH LABEL

tubo para orina con etiqueta / tube pour urine avec étiquette / provetta per urina con etichetta

STERILE

estéril / stérile / sterile



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYSTYRENE

poliestireno / polystyrène / polistirolo

CODE	CAPACITY	PACKAGING
25190	12 ml	INDIVIDUALLY WRAPPED, IN SINGLE BOX en bolsa individual, en estuche / emballage individuel, en étui / in sacchetto singolo, in astuccio

24 H URINE TANK

tanque para orina de 24 h / tank pour urine de 24 h / tanica per raccolta urina 24 h

GRADUATED

graduado / gradué / graduato

WITH SCREW CAP

con tapón roscado / avec bouchon à vis jointé / con tappo a vite



POLYPROPYLENE

polipropileno / polypropylène / polipropilene



POLYETHYLENE

polietileno / polyéthylène / polietilene

CODE	CAPACITY	PACKAGING
25106	2500 ml	IN SINGLE BOX en estuche / en étui / in astuccio

