

## VACUMED® GEL SEPARATOR + LITHIUM HEPARIN

Vacumed® gel separador + Heparina de Litio / Vacumed® gel de séparation + Héparine de Lithium / Vacumed® gel separatore + Litio Eparina

**TORQUOISE CAP**  
tapón turquesa /  
bouchon turquoise /  
tappo turchese

**STERILE**  
estéril /  
stérile /  
sterile

**PET TUBES**  
tubos en PET /  
tubes en PET /  
provette in PET



**42810**



13x75 mm

**2 ml**  
of blood / de sangre /  
de sang / di sangue

**42811**



13x75 mm

**3 ml**  
of blood / de sangre /  
de sang / di sangue

**43816**



13x100 mm

**5 ml**  
of blood / de sangre /  
de sang / di sangue

**44818**



16x100 mm

**8 ml**  
of blood / de sangre /  
de sang / di sangue

## VACUMED® PLUS GEL SEPARATOR + LITHIUM HEPARIN

Vacumed® Plus gel separador + Heparina de Litio / Vacumed® Plus gel de séparation + Héparine de Lithium / Vacumed® Plus gel separatore + Litio Eparina

**TORQUOISE CAP**  
tapón turquesa /  
bouchon turquoise /  
tappo turchese

**STERILE**  
estéril /  
stérile /  
sterile

**PET TUBES**  
tubos en PET /  
tubes en PET /  
provette in PET

**TRANSPARENT LABEL**  
etiqueta transparente /  
étiquette transparente /  
etichetta trasparente



**46810**



13x75 mm

**2 ml**  
of blood / de sangre /  
de sang / di sangue

**46811**



13x75 mm

**3 ml**  
of blood / de sangre /  
de sang / di sangue

**47816**



13x100 mm

**5 ml**  
of blood / de sangre /  
de sang / di sangue

**48818**



16x100 mm

**8 ml**  
of blood / de sangre /  
de sang / di sangue