

## VACUMED® FOR ESR

Vacumed® para VSG / Vacumed® pour VS / Vacumed® per VES

**BLACK CAP**

tapón negro /  
bouchon noir /  
tappo nero

**STERILE**

estéril /  
stérile /  
sterile

**PET TUBES**

tubos en PET /  
tubes en PET /  
provette in PET



**42508**



13x75 mm

**3,2 % x 1,6 ml**

of blood / de sangre /  
de sang / di sangue

**42510**



13x75 mm

**3,2 % x 1,8 ml**

of blood / de sangre /  
de sang / di sangue

**42512**



13x75 mm

**3,2 % x 2,4 ml**

of blood / de sangre /  
de sang / di sangue

## VACUMED® VIRUS STABILIZATION

Vacumed® estabilización de virus / Vacumed® stabilisation du virus /  
Vacumed® stabilizzazione virus

**RED CAP**

tapón rojo /  
bouchon rouge /  
tappo rosso

**STERILE**

estéril /  
stérile /  
sterile

**PET TUBES**

tubos en PET /  
tubes en PET /  
provette in PET



**43500**

**2 ml**

of PBS / de PBS /  
de PBS / di PBS

13x100 mm