

URINTRANSFER® URINE KITS

Urintransfer® kits para orina / Urintransfer® kits pour urine / Urintransfer® kit urina

- YELLOW CAP / tapón amarillo / bouchon jaune / tappo giallo
- STERILE / estéril / stérile / sterile
- PAPER LABEL / etiqueta de papel / étiquette en papier / etichetta cartacea
- POLYPROPYLENE CONTAINER / frasco en polipropileno / pot en polypropylène / contenitore in polipropilene

- PET TUBES / tubos en PET / tubes en PET / provette in PET
- POLYETHYLENE CAP / tapón en polietileno / bouchon en polyéthylène / tappo in polietilene



45080



60 ml (container)
10 ml (test tube)

KIT URINE CONTAINER + VACUUM TUBE

kit frasco para orina + tubo de vacío / kit pot pour urine + tube sous vide / kit contenitore per urina + provetta sottovuoto

45070



120 ml (container)
10 ml (test tube)

KIT URINE CONTAINER + VACUUM TUBE

kit frasco para orina + tubo de vacío / kit pot pour urine + tube sous vide / kit contenitore per urina + provetta sottovuoto

25178



120 ml

URINE CONTAINER FOR VACUUM TUBES

frasco para orina para sistema de vacío / pot pour urine pour système sous vide / contenitore per urina con sistema di aspirazione per provette sottovuoto

25198



10 ml

VACUUM TUBE WITHOUT PRESERVATIVE

tubo de vacío sin preservativo / tube sous vide sans preservatif / provetta sottovuoto senza preservante

25179



120 ml (container)
10 ml (test tube)

KIT URINE CONTAINER + VACUUM TUBE

kit frasco para orina + tubo de vacío / kit pot pour urine + tube sous vide / kit contenitore per urina + provetta sottovuoto