


SEDI-RATE SYSTEM

Sistema Sedi-Rate / Système Sedi-Rate / Sistema Sedi-Rate

TAKIVES SYSTEM

Sistema Takives / Système Takives / Sistema Takives

50334



SEDI-RATE PIPETTE

pipeta Sedi-Rate /
pipette Sedi-Rate /
pipetta Sedi-Rate

GRADUATED FROM 0 TO 180 MM

graduada de 0 a 180 mm /
gradué de 0 à 180 mm /
graduata da 0 a 180 mm

POLYSTYRENE

poliestireno /
polystyrène /
polistirolo

50315



TAKIVES IN GLASS PIPETTE WITH AUTOMATIC ASPIRATOR

pipeta Takives en
vidrio con aspirador
automático / pipette
Takives en verre avec
aspirateur automatique /
pipetta Takives in
vetro con aspiratore
automatico

GRADUATED FROM 0 TO 180 MM

graduada de 0 a 180 mm /
gradué de 0 à 180 mm /
graduata da 0 a 180 mm

GLASS

vidrio /
verre /
vetro

50336




SEDI-RATE RACK WITH 10 NUMBERED PLACES

Sedi-Rate gradilla con 10
puestos numerados /
Sedi-Rate portoir à 10
places numérotées /
supporto Sedi-Rate a 10
posti numerati

POLYSTYRENE

poliestireno /
polystyrène /
polistirolo

50320



TAKIVES RACK WITH 10 NUMBERED PLACES

Takives gradilla con 10
puestos numerados /
Takives portoir à 10
places numérotées /
supporto Takives a 10
posti numerati

ABS

ABS / ABS / ABS